


## YASARGIL VESSEL CLIP; STANDARD

permanent; 17mm; opening 11.9mm

180 gms: angled



**SKU:** 244-410-17

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.47

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.47

**SPECIFICATIONS** [angled](#), [permanent](#)