


YASARGIL VESSEL CLIP; STANDARD

permanent; 20mm; opening 13mm

180 gms: angled



SKU: 244-410-20

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE798

AUTHORS [YASARGIL](#)

Medicon 58.55.46

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE798

AUTHORS [YASARGIL](#)

Medicon 58.55.46

SPECIFICATIONS

[angled](#), [permanent](#)