


## YASARGIL VESSEL CLIP; STANDARD

permanent; 20mm; opening 13mm

180 gms: angled



**SKU:** 244-410-20

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE798

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.46

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE798

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.46

## SPECIFICATIONS

[angled](#), [permanent](#)