


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 7.8mm

180 gms; straight



SKU: 244-400-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE760

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE760

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)