


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 5.6mm

200 gms; angled



SKU: 244-413-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE820

AUTHORS [YASARGIL](#)

Medicon 58.55.49

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE820

AUTHORS [YASARGIL](#)

Medicon 58.55.49

SPECIFICATIONS

[angled](#), [permanent](#)