


## YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 6.2mm

200 gms; curved



**SKU:** 244-423-09

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.42

**SPECIFICATIONS** [curved](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.42

**SPECIFICATIONS** [curved](#), [permanent](#)