


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 8.4mm

150 gms; straight; fenestrated Ø 3.5mm




SKU: 244-430-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FE610 |
| AUTHORS | YASARGIL |
| Medicon | 58.53.51 |
| SPECIFICATIONS | permanent , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FE610 |
| AUTHORS | YASARGIL |
| Medicon | 58.53.51 |

SPECIFICATIONS

[permanent](#), [straight](#)