


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 8.4mm

150 gms; straight; fenestrated Ø 3.5mm



SKU: 244-430-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

Aesculap FE610

AUTHORS [YASARGIL](#)

Medicon 58.53.51

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE610

AUTHORS [YASARGIL](#)

Medicon 58.53.51

SPECIFICATIONS

[permanent](#), [straight](#)