


YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 4.9mm

150 gms; angled; fenestrated \varnothing 3.5mm



SKU: 244-433-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.56

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.56

SPECIFICATIONS [angled](#), [permanent](#)