


YASARGIL VESSEL CLIP; STANDARD

permanent; 17mm; opening 10.6mm

200 gms; straight



SKU: 244-400-17

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | permanent , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | permanent , straight |