


## YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.9mm

150 gms; straight; fenestrated Ø 5mm




**SKU:** 244-436-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE640
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.53.60
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>

**Schließkraft**  
**Closing Force**



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE640
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.53.60

## SPECIFICATIONS

[permanent](#), [straight](#)