


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 9.1mm

150 gms; straight; fenestrated Ø 5mm




SKU: 244-436-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE650
AUTHORS	YASARGIL
Medicon	58.53.61
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE650
AUTHORS	YASARGIL
Medicon	58.53.61

SPECIFICATIONS

[permanent](#), [straight](#)