


## YASARGIL VESSEL CLIP; STANDARD

permanent; 12mm; opening 10.3mm

180 gms; straight; fenestrated Ø 5mm



**SKU:** 244-436-12

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**Aesculap** FE660

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.62

**SPECIFICATIONS** [permanent](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE660

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.62

## SPECIFICATIONS

[permanent](#), [straight](#)