


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 7.8mm

180 gms: angled; fenestrated Ø 5mm



SKU: 244-437-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.65

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.65

SPECIFICATIONS [angled](#), [permanent](#)