


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 5.5mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-439-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.67

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.67

SPECIFICATIONS [angled](#), [permanent](#)