


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 5.5mm

180 gms: angled; fenestrated Ø 5mm



SKU: 244-439-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.68

SPECIFICATIONS [angled](#), [permanent](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.68

SPECIFICATIONS [angled](#), [permanent](#)