


YASARGIL VESSEL CLIP; STANDARD

permanent; 25mm; opening 13.3mm

200 gms; straight



SKU: 244-400-25

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE840
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE840
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight