


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.8mm

150 gms; straight; fenestrated Ø 5mm



SKU: 244-445-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.82

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.82

SPECIFICATIONS [permanent](#), [straight](#)