


## YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.8mm

150 gms; straight; fenestrated Ø 5mm



**SKU:** 244-445-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.82

**SPECIFICATIONS** [permanent](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.82

**SPECIFICATIONS** [permanent](#), [straight](#)