


## YASARGIL VESSEL CLIP; STANDARD

permanent; 10mm; opening 8.8mm

150 gms; straight; fenestrated Ø 5mm



**SKU:** 244-445-10

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.83

**SPECIFICATIONS** [permanent](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.83

**SPECIFICATIONS** [permanent](#), [straight](#)