


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 6mm

150 gms; curved



SKU: 244-401-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE742

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE742

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)