


YASARGIL VESSEL CLIP; STANDARD

temporary; 9mm; opening 7mm

90 gms; straight



SKU: 244-500-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE751
AUTHORS	YASARGIL
Medicon	58.54.60
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE751
AUTHORS	YASARGIL
Medicon	58.54.60

SPECIFICATIONS

[straight](#), [temporary](#)