


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 6.8mm

180 gms; curved



SKU: 244-401-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE752

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE752

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)