


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 5.5mm

110 gms; curved



SKU: 244-503-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.72

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.72

SPECIFICATIONS [curved](#), [temporary](#)