


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 8mm

110 gms; angled



SKU: 244-510-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT270T
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT270T
MATERIAL	Titanium
AUTHORS	YASARGIL

SPECIFICATIONS

[angled](#), [temporary](#)