


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 20mm; opening 13mm

90 gms; angled



SKU: 244-510-20T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.63.20
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.63.20

SPECIFICATIONS

[angled](#), [temporary](#)