


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 8mm; opening 5.4mm

130 gms; curved




SKU: 244-515-08T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.08
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.08

SPECIFICATIONS

[curved](#), [temporary](#)