


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 5.7mm

90 gms; straight




SKU: 244-521-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.37
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.37

SPECIFICATIONS

[straight](#), [temporary](#)