


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 6.2mm

90 gms; straight



SKU: 244-523-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.49
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.65.49

SPECIFICATIONS

[straight](#), [temporary](#)