


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 12mm; opening 9.6mm

90 gms; straight; fenestrated Ø 3.5mm



SKU: 244-530-12T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.12

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.12

SPECIFICATIONS

[straight](#), [temporary](#)