


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

150 gms; curved



SKU: 244-403-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT744T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.28
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT744T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.61.28

SPECIFICATIONS

[curved](#), [permanent](#)