


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-541-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.05

SPECIFICATIONS [angled](#), [temporary](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.05

SPECIFICATIONS

[angled](#), [temporary](#)