


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 6.5mm

180 gms; curved




SKU: 244-403-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.32 |
| SPECIFICATIONS | curved , permanent |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.32 |

SPECIFICATIONS

[curved](#), [permanent](#)