


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 5.6mm

110 gms; angled; fenestrated Ø 5mm



**SKU:** 244-553-05T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.47.65

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.47.65

## SPECIFICATIONS

[angled](#), [temporary](#)