


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 5.6mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-553-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.67

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.67

SPECIFICATIONS

[angled](#), [temporary](#)