


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 15mm; opening 7.4mm

200 gms; curved



SKU: 244-403-15T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.36
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.36

SPECIFICATIONS

[curved](#), [permanent](#)