


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 5.4mm

200 gms; curved




SKU: 244-405-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.37
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.37

SPECIFICATIONS

[curved](#), [permanent](#)