


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 17mm; opening 11.9mm

180 gms: angled



SKU: 244-410-17T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.47

SPECIFICATIONS [angled](#), [permanent](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.47

SPECIFICATIONS

[angled](#), [permanent](#)