


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 6.2mm

200 gms; angled



SKU: 244-411-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


Aesculap FT603T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.50

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT603T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.50

SPECIFICATIONS

[angled](#), [permanent](#)