


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.5mm

200 gms; angled



**SKU:** 244-411-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT613T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.61.51

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT613T

**MATERIAL** [Titanium](#)

## AUTHORS

[YASARGIL](#)

## Medicon

58.61.51

## SPECIFICATIONS

[angled](#), [permanent](#)