


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.5mm

200 gms; angled




SKU: 244-411-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT613T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.51
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT613T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.61.51

SPECIFICATIONS

[angled](#), [permanent](#)