


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 7.8mm

200 gms; angled




SKU: 244-411-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.48
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.48

SPECIFICATIONS

[angled](#), [permanent](#)