


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 7.8mm

180 gms; straight



**SKU:** 244-400-11T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT760T
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.60.12
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT760T
<b>MATERIAL</b>	<a href="#">Titanium</a>

## AUTHORS

[YASARGIL](#)

## Medicon

58.60.12

## SPECIFICATIONS

[permanent](#), [straight](#)