


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 7.8mm

180 gms; straight



SKU: 244-400-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT760T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.12
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT760T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.60.12

SPECIFICATIONS

[permanent](#), [straight](#)