


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 7.2mm

200 gms; curved



SKU: 244-417-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.65
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.65

SPECIFICATIONS

[curved](#), [permanent](#)