


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 8.2mm

200 gms; curved



SKU: 244-419-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT833T


MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.67

SPECIFICATIONS [curved](#), [permanent](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT833T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.67

SPECIFICATIONS

[curved](#), [permanent](#)