


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 6.2mm

200 gms; curved



**SKU:** 244-423-09T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT747T
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.61.42
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT747T
<b>MATERIAL</b>	<a href="#">Titanium</a>

**AUTHORS**

[YASARGIL](#)

**Medicon**

58.61.42

**SPECIFICATIONS**

[curved](#), [permanent](#)