


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 13mm; opening 7.5mm

200 gms; angled




SKU: 244-425-13T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.96 |
| SPECIFICATIONS | angled , permanent |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.96 |

SPECIFICATIONS

[angled](#), [permanent](#)