


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 5.8mm

150 gms; angled; fenestrated Ø 3.5mm



SKU: 244-431-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.53
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.53

SPECIFICATIONS

[angled](#), [permanent](#)