


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 7.2mm

180 gms: angled; fenestrated Ø 3.5mm



SKU: 244-431-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.55
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.55

SPECIFICATIONS

[angled](#), [permanent](#)